**EN** DeVilbiss DV6HH Series CPAP  
Humidifier Use and Care  
Instructions

RX ONLY

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Made in USA of US and Imported Parts.

## IMPORTANT SAFEGUARDS

- Refer to International Standard IEC 60601-1 Ed 3.0 Amendment 1 for safety requirements applicable to Medical Electrical Systems.
- The height of the DV6 series device must be lower than the mask when using a humidifier to prevent water from getting into the mask.
- The humidifier water chamber is intended for single-patient use only.
- Empty and dry humidifier water chamber before transporting.





- Do not attempt to fill the water chamber while it is attached to the humidifier cradle. Damage to the humidifier cradle may occur. Always remove the water chamber from the humidifier cradle before filling.
- Never touch the heater plate on the humidifier cradle. Never touch the heat transfer plate on the bottom of the water chamber. These plates can reach temperatures as high as 158°F (70°C) during operation.
- Do not operate the heater if the water chamber is empty. The heater plate may be turned off using CPAP Settings when device is being used without water.

## INTRODUCTION

### Indications For Use - Optional Heated Humidifier

Use on the advice and prescription of a licensed physician to help relieve the symptoms of dryness of the throat, nasal passages and the mouth, which are common with positive airway pressure therapy. This is especially true in dry climates and during the cold season when humidity in the air is typically lower than at other times.

## SYMBOL DEFINITIONS

	Consult instructions for use		RTCA / DO-160
	Rx Only = Prescription required		TUV Rheinland C-US approval mark
	Class II electrical protection-double insulated		Catalog Number
	Type BF equipment-applied part		Serial Number
	Manufacturer		European Representative
	Date of Manufacture		European Rep CE mark
	Ingress Protection - Protected against finger access to hazardous parts; protected against vertically falling water drops		
	This device contains electrical and/or electronic equipment that must be recycled per EU Directive 2012/19/EU-Waste Electrical & Electronic Equipment		

## IMPORTANT PARTS

### Contents

All models with humidification include one cradle option:

OR

Heated Humidifier Cradle with Standard Humidification (#DV6HH)

Heated Humidifier Cradle with PulseDose® Humidification (#DV6HHPD)

Humidifier Water Chamber (#DV6C)


### Other Replacement Items


Chamber Sealing Gasket (#DV6C-614)

## Understanding your Device's Humidification Options (if prescribed)

### Humidifier Cradle

### Water Chamber

## SYSTEM ASSEMBLY


1. Check to be sure Air Supply Port Cover is inserted into back of CPAP device.

2. Place front of CPAP device on humidifier cradle, ensure locking tabs are engaged, and click into place.

3. Fill clean chamber to fill line with distilled water and insert into cradle. **NOTE** – Use distilled water to avoid mineral buildup inside the chamber.

Chamber Fill


Chamber in Cradle

4. Insert power cord into back of device and wall outlet.

5. Attach tubing to device.

6. Prepare mask and connect to tubing.


## Humidity Settings

1. Select Humidity

2. Choose humidity setting of OFF or 1 through 5 bars. **NOTE:** 5 bars is the highest humidity.


## CLEANING

### Humidifier Cradle

The humidifier cradle should be cleaned weekly.

1. Turn off the heater and unplug from power source. Allow system to cool for at least 10 minutes. Disconnect tubing.

2. Detach device from cradle; first remove water chamber then press flow generator release button. **NOTE** – Flow generator release button must be pressed completely to separate.

3. Wipe the cradle with clean damp cloth.

### Humidifier Water Chamber

Humidifier water chamber should be cleaned daily.


1. Press chamber release latch and remove chamber from cradle.
4. Allow the parts to dry before reassembly.

2. Press forward on chamber disassembly lever and lift the lid to separate halves. **NOTE** – Sealing gasket does not need to be removed for cleaning.

5. Place the tabs on the chamber lid into the slots at the rear of the chamber base and press front until it clicks into place.

3. Hand wash in a solution of warm water and mild dishwashing detergent (e.g. Dawn® dishwashing liquid).


## MAINTENANCE

### Sealing Gasket Placement

**CAUTION** – Inspect the sealing gasket part # DV6C-614 before each use to ensure that it is properly seated under the spokes and in the groove around the entire chamber.

1. Reattach sealing gasket if detached or removed for cleaning. Ensure it snaps under the spokes.

**CAUTION** – Do not use solutions that contain bleach, chlorine, alcohol, scented oils or other aromatic-based ingredients to clean the humidifier components or tubing. These solutions may reduce the life of the product.

**NOTE** – In addition to hand washing, the disassembled chamber may also be washed once per week in a dishwasher, top shelf only.

## MAINTENANCE

**Humidifier Water Chamber** (single-patient use) - Clean the chamber daily. The chamber should be replaced according to the schedule recommended by your provider. Inspect the chamber monthly; replace if damaged.

**Humidifier Cradle** - No maintenance needed for humidifier cradle.

## EXPECTED SERVICE LIFE

Humidifier Cradle - 5 years

Humidifier PulseDose Valve (only model DV6HHPD) - 5 years

Humidifier Water Chamber - 1 year

## TROUBLESHOOTING

Issue	Possible Cause	Remedy
<b>Water condensation collecting in the hose causing a gurgling sound when using optional humidifier.</b>	1. Humidification is too high.	1. Reduce humidifier heater setting.
	2. Room temperature fluctuates from higher to lower levels through the night.	2. Increase room temperature.
<b>Heater setting is not visible on display.</b>	1. The flow generator is not fully seated onto the humidifier cradle.	1. Refer to System Assembly to ensure proper contact is made.

## SPECIFICATIONS

### Humidifier

Size with CPAP Device: .....	6.3" (16.1 cm) H x 6.1" (15.5 cm) W x 7.6" (19.3 cm) D
Weight: .....	3.8 lbs (1.75 kg) w/CPAP device
Humidity Output (in the operating flow range): .....	≥10 mgH <sub>2</sub> O/l air
Maximum Power Consumption: .....	40 watts (humidifier)
Typical Power Consumption with CPAP .....	25 watts
DV6HH Electrical Requirements: .....	11-17 VDC powered from CPAP
Heater Plate Thermal Cutout: .....	239°F (115°C)
Power/Temperature Control: .....	Humidity Setting of 1 (minimum) to 5 (maximum)
Heater Plate Temperature: .....	84°F to 158°F (approx. 29°C to approx. 70°C)
Operating Temperature Range: .....	41°F to 104°F (5°C to 40°C)
Operating Humidity Range: .....	15% to 93% R.H. non-condensing
Operating Atmospheric Pressure Range: .....	70.0 – 106.0 kPa
Operating Flow Range: .....	15 L/M to 60 L/M
Not manufactured with natural rubber latex.	

## PRODUCT DISPOSAL

The outer packaging is made of environmentally friendly materials that can be used as secondary raw materials. If you no longer need this packaging, bring it to your local recycling and waste disposal facility according to the applicable regulations.

The equipment, including accessories and internal components, does not belong in your regular household waste; such equipment is manufactured from high-grade materials and can be recycled and reused. The PC board, heater plate, and wiring should be removed from the unit and recycled as electronic waste. The remaining plastic components should be recycled as plastics.

The European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) requires that electrical and electronic equipment be collected and disposed of separately from other unsorted municipal waste, with the aim of recycling it. The crossed out waste bin symbol indicates that separate collection is required.

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DeVilbiss® and PulseDose® are registered trademarks of DeVilbiss Healthcare. Dawn® is a registered trademark of Procter & Gamble.
© 2015 DeVilbiss Healthcare LLC. 10.15 All Rights Reserved.
**DeVilbiss Healthcare LLC**
**DeVilbiss Healthcare Ltd**
**DeVilbiss Healthcare SAS**
**DeVilbiss Healthcare Pty. Limited**