



**ResMed**

**AirTouch™** **N30i**  
Nasal mask

**AirFit™** **N30i**  
Nasal mask

**AirFit™** **P30i**  
Nasal pillows mask

Frame Fitting Instructions

|                                                                                                       |                                                                                                          |
|-------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------|
| <br>Small<br>Frame | <br>Standard<br>Frame |
| 0–150mm                                                                                               | >150mm                                                                                                   |
| 0–6"                                                                                                  | >6"                                                                                                      |

